

Mouse Nociceptin receptor α OPRL1 α ELISA Kit

Catalog No: #EK20801



Package Size: #EK20801-1 48T #EK20801-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Nociceptin receptor α OPRL1 α ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78-5000pg/mL

Sensitivity:37.5pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.