

Human 39S ribosomal protein L23, mitochondrial (MRPL23) ELISA Kit

Catalog No: #EK20813



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human 39S ribosomal protein L23, mitochondrial (MRPL23) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |
| Other Names | FLJ45387; L23MRP; RPL23; RPL23L; L23 mitochondrial-related protein OTTHUMP00000069610 ribosomal protein L23-like ribosomal protein related to L23 (mitochondrial) |

Application Details

| |
|--|
| Detect Range:Request Information |
| Sensitivity:Request Information |
| Sample Type:Serum, Plasma, Other biological fluids |
| Sample Volume: 1-200 μL |
| Assay Time:1-4.5h |
| Detection wavelength:450 nm |

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.