

# Human TRAF-type zinc finger domain-containing protein 1(TRAFF1) ELISA Kit

Catalog No: #EK20842

Package Size: #EK20842-1 48T #EK20842-2 96T



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human TRAF-type zinc finger domain-containing protein 1(TRAFF1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.069ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.