

Mouse Chemokine-like protein TAFA-4(TAFA4) ELISA Kit

Catalog No: #EK20818

Package Size: #EK20818-1 48T #EK20818-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Chemokine-like protein TAFA-4(TAFA4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:12.5pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.