

ζ δ1 ι ε 1ε ζA²(2%)

Catalog No: #AR0735



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ζ δ1 ι ε 1ε ζA ² (2%)
Storage	4°C, 1.0 Ωε , 12δE ζ

Application Details

ε 1ε η η η ε "δ1 δΩ• ε- η» θ ε η» η» η ε ζ ε1/2ζ • η» ζ οΟ ε 1ε ε θ 1/2 »ζ-ε ζ1• ζ §ζ1Α1 Αδ1/2 ι Αε—Hθ θ1«η
 » η» ε η» θ η θ ζ1Αγ • ζ ε Αη» θ• ε ι θ• η η ι ΩοΟ ε 1ε ε ι θΩ ε • ε 1γ • η ζ • ζ»ε ε ε• η ©η- δ1/2Ωθ
 η 1/2θ"ε ι "η» ζ ε• η ζ Ηε• οΟ δ» θ δ1/2Ωι Αε ±ζ'»γ ε 1ε ζA²δE»θ1• ε δE1ι η±»ε 1ε ζA²γ • ζ± η±»ε 1ε ζA²γ • ι η
 ±»ε 1ε ζA²γ • ζ°§ε ε η±»ε 1ε ζA²γ • θ 1ε ρι 'Eη • η±»ε 1ε ζA²η- οΟ θY δE1ε'EΕη "η ζ —ι η±»δE-η η1• ε° ι©-ζ γ • ι η±
 »δE-η δH ι γ ζ δ1 ι ε 1ε ζA²δO εO θ"θ η 1/2θ"α-θ 1ζ η» ζ ε• ε1/2οO δE• ε ©δ1 θΩ ζ°§ε η ©ι Αζ θ 2οO ι ι ε1'
 εΩ«οO ζE ι • ε ε• γ

ζ δ1 ι ε 1ε ζA²(2%)δE»θ1• η ±ζ δ1 ι γ • η£-ι 'Eη • η- η» ζ • οO pHε 'O7.2ο1/2 7.4οO θ— ε 1ε ζA²ε—Hη» θ ζ'Eγ • η» θ ζ'
 η η» εY η» ζ ε 1ε ζ ζ ε 1/2οO ζ —ζ ε'EΕη "η ζ ε ζ δ1 ι ε 1ε ζA²οO η»• ε'EΕη "δ1 η "ι ζ ζ η-η ε 1ε γ θ
 — θ— ε δ» η "δ1 η§ η ι' ε οO δE• ι η "δ1 δE'ε1 θ— ζ -ζ ε Αδ» η "ι γ

Note: This product is for in vitro research use only and is not intended for use in humans or animals.