

# PANK3 rabbit pAb

Catalog No: #62057



Package Size: #62057-1 50ul #62057-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PANK3 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human Mouse
Specificity	This antibody detects endogenous levels of human PANK3
Immunogen Description	Synthesized peptide derived from human PANK3
Calculated MW	41kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.