

Junctophilin 4 Antibody HRP Conjugated

Catalog No: #C00936H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Junctophilin 4 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Junctophilin 4
Conjugates	HRP
Target Name	Junctophilin 4
Other Names	JP-4; Jph4; JPH4_HUMAN; JPHL1; Junctophilin4; Junctophilin like 1 protein; Junctophilin-4; Junctophilin-like 1 protein; KIAA1831.
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Junctophilin 4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.