

## JAMC Antibody HRP Conjugated

Catalog No: #C00938H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JAMC Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 45-95 310 derived from human JAMC
Conjugates	HRP
Target Name	JAMC
Other Names	JAMC; JAM-2; JAM-3; JAM-C; Junctional adhesion molecule C; Junctional adhesion molecule 3; JAM3; UNQ859 PRO1868
Accession No.	Swiss-Prot#Q9BX67NCBI Gene ID83700
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

## Background

Participates in cell-cell adhesion. It is a counter-receptor for ITGAM, mediating leukocyte-platelet interactions and is involved in the regulation of transepithelial migration of polymorphonuclear neutrophils (PMN). The soluble form is a mediator of angiogenesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.