

# Ninjurin 1 Antibody Biotin Conjugated

Catalog No: #C00953B

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ninjurin 1 Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-70 152 derived from human Ninjurin 1
Conjugates	Biotin
Target Name	Ninjurin 1
Other Names	NIN1; NINJURIN; Ninjurin-1; Nerve injury-induced protein 1; NINJ1
Accession No.	Swiss-Prot#Q92982NCBI Gene ID4814
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000

## Background

Homophilic cell adhesion molecule that promotes axonal growth. May play a role in nerve regeneration and in the formation and function of other tissues. Cell adhesion requires divalent cations.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.