

JunB (Phospho-Thr102 Thr104) Antibody Biotin Conjugated

Catalog No: #C01031B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JunB (Phospho-Thr102 Thr104) Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human JunB around the phosphorylation site of Thr102 Thr104 [TT(p-T)P(p-T)PP]
Conjugates	Biotin
Target Name	JunB Thr102 + Thr104
Other Names	Transcription factor jun-B; JUNB
Accession No.	NCBI Gene ID3726
Concentration	1mg/ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.