

WFS1 Antibody HRP Conjugated

Catalog No: #C01071H

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	WFS1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 780-830 890 derived from human WFS1
Conjugates	HRP
Target Name	WFS1
Other Names	WFS; WFRS; WFSL; CTRCT41; Wolframin; WFS1
Accession No.	Swiss-Prot#O76024NCBI Gene ID7466
Cell Localization	Lumenal
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Participates in the regulation of cellular Ca(2+) homeostasis, at least partly, by modulating the filling state of the endoplasmic reticulum Ca(2+) store.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.