

VKORC1 Antibody


Catalog No: #25474

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VKORC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the amino terminus of human VKORC1.
Target Name	VKORC1
Other Names	Vitamin K epoxide reductase complex subunit 1, VKOR, MST134, MST576, VKCFD2, EDTP308
Accession No.	NP_076869
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of VKORC1 in A549 cell lysate with VKORC1 antibody at 1 ug/mL.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.