

# TAB1(N-term) Monoclonal Antibody

Catalog No: #27075

Package Size: #27075-1 50ul #27075-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	TAB1(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu D Ms P Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human TAB1
Target Name	TAB1
Other Names	TGF-beta-activated kinase 1 and MAP3K7-binding protein 1 □ TAK1-binding protein 1 □ transforming growth factor beta-activated kinase-binding protein 1 □ mitogen-activated protein kinase kinase kinase 7-interacting protein 1
Accession No.	UniProt:Q15750 NCBI Gene ID:10454
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 55kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected TAB1 on 40ug of cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.