

## PPP5C Monoclonal Antibody

Catalog No: #27096

Package Size: #27096-1 50ul #27096-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PPP5C Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2F2-F7-E8
Isotype	IgG2a
Applications	ELISA WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human PPP5C protein fragments expressed in E.coli.
Target Name	PPP5C
Other Names	PP5, PPT, PPP5
Formulation	ascites

## Application Details

Predicted MW: 57kd

Western blotting: 1:1000

## Images



Western blot detection of PPP5C in HeLa

## Background

May play a role in the regulation of RNA biogenesis and/or mitosis. In vitro, dephosphorylates serine residues of skeletal muscle phosphorylase and histone H1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.