

NUP98 Monoclonal Antibody

Catalog No: #27111



Package Size: #27111-1 50ul #27111-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	NUP98 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B8-D7-H10
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Target Name	NUP98
Other Names	ADIR2, NUP96, NUP196
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Application Details


Predicted MW: 98kd

Western blotting: 1:1000

Images



Western blot detection of NUP98 in HeLa



Immunoprecipitation analysis of HeLa cell lysates using NUP98 antibody

Background

Nup98 and Nup96 play a role in the bidirectional transport across the nucleoporin complex (NPC). The repeat domain in Nup98 has a direct role in the transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.