

JHDM1 ANTIBODY HRP Conjugated

Catalog No: #C02711H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JHDM1 ANTIBODY HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from <i>Saccharomyces cerevisiae</i> YJM789 of JHDM1
Conjugates	HRP
Target Name	JHDM1
Other Names	NULL
Accession No.	NCBI Gene ID856777
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

KDM2 JHD1 is a JmjC domain family histone demethylase specific for H3-K36, similar to proteins found in human, mouse, drosophila, *X. laevis*, *C. elegans*, and *S. pombe*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.