

WASP Antibody HRP Conjugated

Catalog No: #C02773H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | WASP Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 210-260 502 derived from human WASP |
| Conjugates | HRP |
| Target Name | WASP |
| Other Names | THC; IMD2; SCNX; THC1; WASP; WASPA; Wiskott-Aldrich syndrome protein; WAS |
| Accession No. | Swiss-Prot#P42768NCBI Gene ID7454 |
| Cell Localization | Cytoplasm |
| Concentration | 1mg ml |
| Formulation | 10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin. |
| Storage | Store at 4C for 12 months. |

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Effector protein for Rho-type GTPases. Regulates actin filament reorganization via its interaction with the Arp2 3 complex. Important for efficient actin polymerization. Possible regulator of lymphocyte and platelet function. Mediates actin filament reorganization and the formation of actin pedestals upon infection by pathogenic bacteria.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.