

WDR91 Antibody Biotin Conjugated

Catalog No: #C02775B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WDR91 Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human WDR91
Conjugates	Biotin
Target Name	WDR91
Other Names	HSPC049; WD repeat containing protein 91; WD repeat domain 91; WD repeat-containing protein 91; WDR91; WDR91; WDR91_HUMAN.
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

WDR91

Note: This product is for in vitro research use only and is not intended for use in humans or animals.