

WNT2B Antibody HRP Conjugated

Catalog No: #C03819H

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	WNT2B Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 341-391 391 derived from human WNT2B
Conjugates	HRP
Target Name	WNT2B
Other Names	WNT13; Protein Wnt-2b; Protein Wnt-13; WNT2B
Accession No.	Swiss-Prot#Q93097NCBI Gene ID7482
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters. May be involved in normal development or differentiation as well as in carcinogenesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.