

## VSV-G tag Antibody HRP Conjugated

Catalog No: #C03953H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VSV-G tag Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	VSV-G tagged protein
Crossing Reactivity	VSV-G tagged protein
Immunogen Description	KLH conjugated VSV-G tag
Conjugates	HRP
Target Name	VSV-G tag
Other Names	Vesicular stomatitis virus glycoprotein tag; VSV epitope tag; VSV tag; VSV-G epitope tag; VSVG; vsvg tag.
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:10000-20000

## Background

Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.