

Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody

Catalog No: #11584




Package Size: #11584-1 50ul #11584-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Applications | WB |
| Species Reactivity | Hu |
| Immunogen Type | Peptide-KLH |
| Target Name | Acetyl-CoA Carboxylase |
| Modification | Phospho |
| Other Names | ACC, ACAC, ACC1, ACCA, |
| Concentration | 1.0mg/ml |

Images



Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody #11584.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.