

MBP-Tag Mouse Monoclonal Antibody

Catalog No: #T601



Package Size: #T601-1 50ul #T601-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MBP-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB IHC IF IP
Specificity	MBP-tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal MBP -tagged recombinant proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Full-length MBP protein
Target Name	MBP-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

Images



Western blot analysis of MBP-tag fusion protein, using MBP-tag mouse monoclonal antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.