

## T7-Tag Mouse Monoclonal Antibody

Catalog No: #T604



Package Size: #T604-1 50ul #T604-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                       |  |
|-----------------------|--|
| Product Name          | T7-Tag Mouse Monoclonal Antibody   |
| Host Species          | Mouse  |
| Clonality             | Monoclonal   |
| Isotype               | IgG1   |
| Purification          | Affinity-chromatography  |
| Applications          | WB IF/ICC  |
| Specificity           | T7-tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal T7-tagged proteins. |
| Immunogen Type        | Peptide  |
| Immunogen Description | A synthetic peptide MASMTGGQQMG  |
| Target Name           | T7-Tag   |
| SDS-PAGE MW           | N/A  |
| Concentration         | 1.0mg/ml   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.  |
| Storage               | Store at -20°C. Stable for 12 months from date of receipt  |

## Application Details

WB: 1:3000~1:10000

IF/ICC: 1:200

## Images



Western blot analysis of GST-T7 fusion protein, using T7-tag mouse monoclonal antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.