

Jak1 (Phospho-Tyr1034 + Tyr1035) Antibody FITC Conjugated

Catalog No: #C04568F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Jak1 (Phospho-Tyr1034 + Tyr1035) Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB RtB B B
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human Jak1 around the phosphorylation site of Tyr1034 + Tyr1035
Conjugates	FITC
Target Name	Jak1 Tyr1034 + Tyr1035
Other Names	Tyrosine-protein kinase JAK1; Janus kinase 1; JAK-1; JAK1; JAK1A; JAK1B
Accession No.	NCBI Gene ID3716
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200B IF=1:50-200

Background

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha beta gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.