

FOXS1 Antibody

Catalog No: #31072



Package Size: #31072-1 50ul #31072-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FOXS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	ELISA WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total FOXS1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	FOXS1
Other Names	forkhead box S1, FKHL18, FREAC10
Accession No.	Genbank No.: BC013408
Concentration	0.8mg/ml
Formulation	Supplied at 0.5mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year


Application Details

Predicted MW: 35kd

ELISA: 1:2000-1:5000

Western blotting: 1:500-1:2000

Images



Gel: 10%+12%SDS-PAGE

Lysate: 40 µg Mouse heart tissue lysate

Primary antibody: 1/250 dilution

Secondary antibody: Goat anti Rabbit IgG - H&L (HRP) at 1/10000 dilution

Exposure time: 10 seconds

Background

The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.