

YAP1(Phospho-Tyr407) Antibody Biotin Conjugated

Catalog No: #C04732B

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | YAP1(Phospho-Tyr407) Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic phosphopeptide derived from human YAP around the phosphorylation site of Tyr407 |
| Conjugates | Biotin |
| Target Name | YAP1 Tyr407 |
| Other Names | 65 kDa Yes associated protein; YAp 1; YAp-1; YAP 65; YAP; YAP2; YAP65; Yes associated protein 1 65kDa; Yes associated protein 1; Yes associated protein 2; YKI. |
| Accession No. | NCBI Gene ID10413 |
| Concentration | 1mg ml |
| Formulation | 10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide. |
| Storage | Store at 4C for 12 months. |

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

YAP1 (also known as Yes associated protein 1) was originally identified as a transcription factor that binds to the SH3 domain of the YES kinase (a Src protein kinase). More recently it has been identified as a candidate oncogene that promotes tumorigenesis in many different types of cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.