

YAP1(Phospho-Tyr407) Antibody HRP Conjugated

Catalog No: #C04732H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	YAP1(Phospho-Tyr407) Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human YAP around the phosphorylation site of Tyr407
Conjugates	HRP
Target Name	YAP1 Tyr407
Other Names	65 kDa Yes associated protein; YAp 1; YAp-1; YAP 65; YAP; YAP2; YAP65; Yes associated protein 1 65kDa; Yes associated protein 1; Yes associated protein 2; YKI.
Accession No.	NCBI Gene ID10413
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

YAP1 (also known as Yes associated protein 1) was originally identified as a transcription factor that binds to the SH3 domain of the YES kinase (a Src protein kinase). More recently it has been identified as a candidate oncogene that promotes tumorigenesis in many different types of cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.