

# WNT10B Antibody Biotin Conjugated

Catalog No: #C04844B

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WNT10B Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 220-270 389 derived from human WNT10b
Conjugates	Biotin
Target Name	WNT10B
Other Names	SHFM6; WNT-12; Protein Wnt-10b; Protein Wnt-12; WNT10B; WNT12
Accession No.	Swiss-Prot#O00744NCBI Gene ID7480
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

## Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters (By similarity).

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.