

WNT10B Antibody HRP Conjugated

Catalog No: #C04844H

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	WNT10B Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 220-270 389 derived from human WNT10b
Conjugates	HRP
Target Name	WNT10B
Other Names	SHFM6; WNT-12; Protein Wnt-10b; Protein Wnt-12; WNT10B; WNT12
Accession No.	Swiss-Prot#O00744NCBI Gene ID7480
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.