

WIF1 Antibody HRP Conjugated

Catalog No: #C05014H

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	WIF1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 379 derived from human WIF1
Conjugates	HRP
Target Name	WIF1
Other Names	WIF-1; Wnt inhibitory factor 1; WIF1; UNQ191 PRO217
Accession No.	Swiss-Prot#Q9Y5W5NCBI Gene ID11197
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

Background

Binds to WNT proteins and inhibits their activities. May be involved in mesoderm segmentation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.