

## JunB(Phospho-Ser259) Antibody Biotin Conjugated

Catalog No: #C05785B

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JunB(Phospho-Ser259) Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human JunB around the phosphorylation site of Ser259 [PV(p-S)PI]
Conjugates	Biotin
Target Name	JunB Ser259
Other Names	Activator protein 1; AP 1; AP1; Jun B; Jun B proto oncogene; Jun B protooncogene; JunB proto oncogene; JunB protoncogene 9; JunB protooncogene; Transcription factor jun B; Transcription factor junB.
Accession No.	NCBI Gene ID3726
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

## Background

The cJun proto-oncogene was first identified as the cellular homolog of the avian sarcoma virus vjun oncogene. JunB and JunD have been shown to be almost identical to cJun in their C terminal regions, which are involved in dimerization and DNA binding, whereas their N terminal domains, which are involved in transcriptional activation, diverge. JunB is a transcription factor involved in regulating gene activity following the primary growth factor response. It binds to the DNA sequence 5'-TGA[CG]TCA-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.