

Wnt8A Antibody Biotin Conjugated

Catalog No: #C06235B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Wnt8A Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 200-250 351 derived from human Wnt8A
Conjugates	Biotin
Target Name	Wnt8A
Other Names	WNT8D; Protein Wnt-8a; Protein Wnt-8d; WNT8A
Accession No.	Swiss-Prot#Q9H1J5NCBI Gene ID7478
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Immunohistochemistry1:100-500

Background

Ligand for members of the frizzled family of seven transmembrane receptors. May play an important role in the development and differentiation of certain forebrain structures, notably the hippocampus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.