

# Wnt8A Antibody HRP Conjugated

Catalog No: #C06235H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Wnt8A Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 200-250 351 derived from human Wnt8A
Conjugates	HRP
Target Name	Wnt8A
Other Names	WNT8D; Protein Wnt-8a; Protein Wnt-8d; WNT8A
Accession No.	Swiss-Prot#Q9H1J5NCBI Gene ID7478
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Immunohistochemistry1:100-500

## Background

Ligand for members of the frizzled family of seven transmembrane receptors. May play an important role in the development and differentiation of certain forebrain structures, notably the hippocampus.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.