

# WNT4 Antibody HRP Conjugated

Catalog No: #C06239H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WNT4 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 351 derived from human WNT4
Conjugates	HRP
Target Name	WNT4
Other Names	WNT-4; SERKAL; Protein Wnt-4; WNT4; UNQ426 PRO864
Accession No.	Swiss-Prot#P56705NCBI Gene ID54361
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

## Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters (By similarity). Overexpression may be associated with abnormal proliferation in human breast tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.