

# WISP39 Antibody Biotin Conjugated

Catalog No: #C06945B

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WISP39 Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human WISP39 FKBPL
Conjugates	Biotin
Target Name	WISP39
Other Names	NG7; WISP39; DIR1; FK506 binding protein like; FK506-binding protein-like; FkbpI; FKBPL_HUMAN; OTTHUMP00000029170; WAF 1 CIP1 stabilizing protein 39; WAF-1 CIP1 stabilizing protein 39.
Accession No.	NCBI Gene ID63943
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500

## Background

May be involved in response to X-ray. Regulates p21 protein stability by binding to Hsp90 and p21. Tissue specificity: Ubiquitously expressed with higher levels in testis.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.