

## TFF3 Antibody

Catalog No: #32532

Package Size: #32532-1 50ul #32532-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TFF3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TFF3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TFF3 (NP_003217.3).
Target Name	TFF3
Other Names	TFF3;ITF;P1B;TFI
Accession No.	Uniprot:Q07654GeneID:7033
SDS-PAGE MW	Refer to figures
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

## Images



Immunofluorescence analysis of human colon carcinoma using TFF3 Rabbit pAb.

## Background

Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa

from insults, stabilize the mucus layer and affect healing of the epithelium. This gene is expressed in goblet cells of the intestines and colon. This gene and two other related trefoil family member genes are found in a cluster on chromosome 21.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.