

## WDR26 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C00502Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WDR26 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P) IF(IHC-F)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 260-310 661 derived from human WDR26
Conjugates	Cy5
Target Name	WDR26
Other Names	CDW2; GID7; MIP2; WD repeat-containing protein 26; CUL4- and DDB1-associated WDR protein 2; Myocardial ischemic preconditioning up-regulated protein 2; WDR26; PRO852
Accession No.	Swiss-Prot#:Q9H7D7NCBI Gene ID:80232
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF:1:50-200

## Background

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Two transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq].

Note: This product is for in vitro research use only and is not intended for use in humans or animals.