

TCN2 Antibody

Catalog No: #33019

Package Size: #33019-1 50ul #33019-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total TCN2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TCN2 (NP_001171655.1).
Target Name	TCN2
Other Names	TCN2;D22S676;D22S750;II;TC;TC II;TC-2;TC2;TCII
Accession No.	Uniprot:P20062GeneID:6948
SDS-PAGE MW	54kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TCN2 antibody.

Background

This gene encodes a member of the vitamin B12-binding protein family. This family of proteins, alternatively referred to as R binders, is expressed in various tissues and secretions. This plasma protein binds cobalamin and mediates the transport of cobalamin into cells. This protein and other

mammalian cobalamin-binding proteins, such as transcobalamin I and gastric intrinsic factor, may have evolved by duplication of a common ancestral gene. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.