

## TOPBP1 Antibody

Catalog No: #33043

Package Size: #33043-1 50ul #33043-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TOPBP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total TOPBP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TOPBP1 (NP_008958.2).
Target Name	TOPBP1
Other Names	TOPBP1;TOP2BP1
Accession No.	Uniprot:Q92547GeneID:11073
SDS-PAGE MW	200KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of HCT116 cells, using TOPBP1 antibody.

## Background

This gene encodes a binding protein which interacts with the C-terminal region of topoisomerase II beta. This interaction suggests a supportive role for this protein in the catalytic reactions of topoisomerase II beta through transient breakages of DNA strands.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.