

# ETV4 Antibody

Catalog No: #33049



Package Size: #33049-1 50ul #33049-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ETV4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total ETV4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ETV4.
Target Name	ETV4
Other Names	E1AF; PEA3; E1A-F; PEAS3;
Accession No.	Swiss-Prot:P43268NCBI Gene ID:2118
SDS-PAGE MW	53KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


## Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using Pea3 / ETV4 antibody at 1:1000 dilution.



Immunofluorescence analysis of U2OS cells using Pea3 / ETV4 antibody. Blue: DAPI for nuclear staining.

## Background

Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.