

## Junctophilin 4 Polyclonal Antibody Cy7 Conjugated

Catalog No: #C00936Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Junctophilin 4 Polyclonal Antibody Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Junctophilin 4
Conjugates	Cy7
Target Name	Junctophilin 4
Other Names	JP-4; Jph4; JPH4_HUMAN; JPHL1; Junctophilin4; Junctophilin like 1 protein; Junctophilin-4; Junctophilin-like 1 protein; KIAA1831.
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF:1:50-200

## Background

Junctophilin 4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.