

# JunB (Phospho-Thr102 Thr104) Polyclonal Antibody Cy5 Conjugated

Catalog No: #C01031Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JunB (Phospho-Thr102 Thr104) Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human JunB around the phosphorylation site of Thr102 Thr104 [TT(p-T)P(p-T)PP]
Conjugates	Cy5
Target Name	JunB Thr102 + Thr104
Other Names	Transcription factor jun-B; JUNB
Accession No.	NCBI Gene ID:3726
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF:1:50-200

## Background

Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.