

## ZNF18 Antibody

Catalog No: #33689

Package Size: #33689-1 50ul #33689-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total ZNF18 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF18.
Target Name	ZNF18
Other Names	Zinc finger protein 18; Zinc finger protein KOX11; Heart development-specific gene 1 protein; Zinc finger protein with KRAB and SCAN domains 6; ZNF18
Accession No.	Swiss-Prot: P17022NCBI Gene ID: 7566
SDS-PAGE MW	62kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from HeLa cells, using ZNF18 antibody #33689.

## Background

---

May be involved in transcriptional regulation.

Guo L.L., Yi Chuan 27:523-530(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Thiesen H.-J., New Biol. 2:363-374(1990).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.