

## ZNF76 Antibody

Catalog No: #33760

Package Size: #33760-1 50ul #33760-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF76 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF76 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF76.
Target Name	ZNF76
Other Names	Zinc finger protein 76; ZNF76; D6S229E;
Accession No.	Swiss-Prot: P36508NCBI Gene ID: 7629
SDS-PAGE MW	62kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:50~1:100

## Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ZNF76 antibody #33760.

## Background

May be involved in transcriptional regulation.

Ragoussis J., Genomics 14:673-679(1992).

Phelan M., Submitted (OCT-2004) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.