

## ZNF134 Antibody

Catalog No: #33807

Package Size: #33807-1 50ul #33807-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZNF134 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF134 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF134.
Target Name	ZNF134
Other Names	Zinc finger protein 134; ZNF134;
Accession No.	Swiss-Prot: P52741NCBI Gene ID: 7693
SDS-PAGE MW	36kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from COS-7 cells and K562 cells, using ZNF134 antibody #33807.

## Background

May be involved in transcriptional regulation.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.