

ZNF596 Antibody

Catalog No: #33986

Package Size: #33986-1 50ul #33986-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ZNF596 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF596 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF596.
Target Name	ZNF596
Other Names	Zinc finger protein 596; ZNF596;
Accession No.	Swiss-Prot: Q8TC21NCBI Gene ID: 169270
SDS-PAGE MW	58kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from Jurkat cells, using ZNF596 antibody #33986.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ZNF596 antibody #33986.



Background

May be involved in transcriptional regulation.

Ota T., Nat. Genet. 36:40-45(2004).

Nusbaum C., Nature 439:331-335(2006).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.