

ZNF541Antibody

Catalog No: #34058

Package Size: #34058-1 50ul #34058-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF541Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF541 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF541.
Target Name	ZNF541
Other Names	Zinc finger protein 541; ZNF541; ZN541;
Accession No.	Swiss-Prot: Q9H0D2NCBI Gene ID: 84215
SDS-PAGE MW	83kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using ZNF541 antibody #34058.

Background

Component of some chromatin remodeling multiprotein complex that plays a role during spermatogenesis By similarity.

Wiemann S., Genome Res. 11:422-435(2001).
Yu L.-R., J. Proteome Res. 6:4150-4162(2007).
Sjoebloom T., Science 314:268-274(2006).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.