

## ZNF95 Antibody

Catalog No: #34129

Package Size: #34129-1 50ul #34129-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF95 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF95 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF95.
Target Name	ZNF95
Other Names	Zinc finger protein with KRAB and SCAN domains 5; Zinc finger protein 95 homolog; Zfp-95; ZKSCAN5; KIAA1015
Accession No.	Swiss-Prot: Q9Y2L8NCBI Gene ID: 23660
SDS-PAGE MW	97kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Immunofluorescence: 1:100~1:500

## Images



Immunofluorescence analysis of A549 cells, using ZNF95 antibody #34129.

## Background

---

May be involved in transcriptional regulation.

Dreyer S.D., Genomics 62:119-122(1999).

Nagase T., DNA Res. 6:63-70(1999).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.