

JHDM1 POLYCLONAL ANTIBODY Cy3 Conjugated

Catalog No: #C02711Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | JHDM1 POLYCLONAL ANTIBODY Cy3 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF(IHC-P) |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide derived from <i>Saccharomyces cerevisiae</i> YJM789 of JHDM1 |
| Conjugates | Cy3 |
| Target Name | JHDM1 |
| Other Names | NULL |
| Accession No. | NCBI Gene ID:856777 |
| Concentration | 1mg ml |
| Formulation | Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide. |
| Storage | Store at 4C for 12 months. |

Application Details

IF:1:50-200

Background

KDM2 JHD1 is a JmjC domain family histone demethylase specific for H3-K36, similar to proteins found in human, mouse, drosophila, *X. laevis*, *C. elegans*, and *S. pombe*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.