

JHDM1 POLYCLONAL ANTIBODY Cy7 Conjugated

Catalog No: #C02711Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JHDM1 POLYCLONAL ANTIBODY Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from <i>Saccharomyces cerevisiae</i> YJM789 of JHDM1
Conjugates	Cy7
Target Name	JHDM1
Other Names	NULL
Accession No.	NCBI Gene ID:856777
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

IF:1:50-200

Background

KDM2 JHD1 is a JmjC domain family histone demethylase specific for H3-K36, similar to proteins found in human, mouse, drosophila, *X. laevis*, *C. elegans*, and *S. pombe*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.